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110°C

400°C

Precision heating baths & Waterbath

for continuous operation

700°C

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PZ26-7

700°C

W16, a water bath for temperatures up to 110°C. The stainless-steel tank (160 mmØ x 100 mm) is tightly mounted in stainless steel housing. A thermostat controls the temperature with an accuracy of ± 1.5 K or better

Dry-boil protection prevents overheating in case of water shortage. A temperature fuse provides additional protection.

The low lying ring-set (Order No. -Nr. RI 17/E) enables the condensed water to drop back into the bath.

PZ26-4, with a bath-capacity of 5 litres, can be used as fluid, air and sandbath for temperatures up to 400°C. The bath filling can be oil, water, sand, or any other filling which will not attack the surface of the bath, which is made out of high-grade, stainless steel (W-No. 1.4301).

PZ26-7, designed for temperatures of up to 700°C offers the same possibilities, but not as a liquid bath. This device can only be used as a dry, air or sand bath. Both baths are also well-suited for rotary evaporators.

PZ26-4 and PZ26-7 are both heated by a tightly soldered tubular heating element.

PZ26-4 und PZ26-7 are microprocessor controlled by a PID Regulator with a low-mass and fast-reacting temperature sensor. The actual temperature is constantly displayed. The set temperature appears at the push of a button and is adjustable with the up/down-button.


The housing is made of high-grade stainless-steel, and the base is additionally painted.

All devices are constructed for continuous operation. Even at maximum temperature, the housing and component parts will heat up only slightly.

W16, P26-4 and PZ26-7 are equipped with a screw-in socket for a support rod 12 mmØ.

Technische Daten Best.-Nr. = Typ		W16	PZ26-4	PZ26-7	Technical Data Order No. = types
Temperaturvorwahl	°C	30...110	20...400	20...700	Range of temperature
Anzeige		ohne/without	digital	digital	Display
Leistung	W	1000	1200	1200	Power
Spannung	V	220...240 V, 50-60 Hz			Voltage
Außenabmessungen	mm	220 x 220	310 x 310		Overall dimensions
Höhe	mm	160	280	280	Height
Gewicht	kg	3,0	6,8	6,8	Weight

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The right to make modification serving further technical development, is reserved

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Further Fabrications

High temperature TITANIUM hotplate up to 600°C - Precision hotplates - Precision hotplate made of cast iron up to 450°C - Stand hotplates with supporting rod - Rapid incinerators - Multipurpose, extremely deep, heating bath - Sandbaths - Series heating units - Multipurpose heating baths - Water- and Oil baths - Development - Special orders.