



Glass door

Allows observation of measured equipment without opening the door.



Cable pass-through

The cable pass-through in the right side wall (7 cm diameter) allows easy connection to devices.

LABO 720 CHROMAT

PRO-ACTIVE





















according to DIN 13277

■ PRO-ACTIVE- Control:

Permanent, proactive monitoring of the performance data and alerts in the case of deviations – so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.

- External housing made from galvanised sheet steel (rustproof), with white, anti-scratch powder coating. Length of the plug cable: approx. 3 m.
- Equipped with castors.
- Interior made from smooth aluminium with colourless protective coating. Gibs (every 15 mm) for flexible interior arrangement.
- Interior consists of five plastic-coated shelves on overlays.
- Extra-thick 70 mm energy saving insulation, made from high-quality, compression-moulded and environmentally friendly material.
- Self-closing glasdoor with easy-to-replace plastic magnetic seal frame, lockable.
- Door hinge on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.
- Forced-air cooling, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviations.
- **LED-Illumination** mounted on side wall. Permanent lighting can be switched-on over a button at the control panel.

- 3 moisture-proof sockets mounted in the interior.
- The cable pass-through in the right side wall (7cm diameter) allows easy connection to devices.
- Automatic fast defrosting thanks to time-limited and temperature-monitored reversal of the refrigerant circuit.
- Condensate evaporates in stainless steel heated dish beneath the refrigerator floor.
- Key switch protected control panel against tampering. Membrane keyboard with digital temperature display and minimum/maximum temperature memory.
- Power failure alarm (visual and audible alarm), the monitoring unit remains in operation for approx. 30 hours on battery power.
- Warning functions with visual and audible alarm signal in the case of temperature deviations or other malfunctions. Door open alarm after 60 seconds.
- Alarms can be forwarded using a potential-free contact (e.g., to a mobile phone with optional GSM-Module or to a control centre).
- Data documentation can be read via USB interface with KIRSCH Datanet software.
- Antifreeze against sub-zero temperatures.
- Ventilation-enforced refrigerating machine, hermetically sealed, energy saving, low-noise, easy to service.
- Low-noise compressor reduces noise to 40 dB (corresponds to noise level I; sound perception: whisper).

Specification

Capacity	700 litres
Temperature setting	approx. +4 °C to +15 °C
Voltage	220 – 240 V, 50/60 Hz
Power consumption	250 watts
Normal consumption	1.63 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	553 watts
Exterior dimensions (including door handle and distance from wall)	w x d x h = 77 x 98 x 195.5 cm
Interior dimensions	w x d x h = 62 x 77 x 140 cm (usable width 2 cm, usable depth 11 cm, usable height 13 cm less)
Exterior dimensions with door open at 90°	w x d = 77 x 165 cm
Shelf size	w x d = 59 x 65 cm
Max. load shelf	40 kg
Weight	Net weight 146 kg, gross weight 174 kg

Optional equipment

- Water-cooled refrigerating machine.
- Adjustable feet to compensate uneven floors.
- Housing made of stainless steel 4301, longitudinal brushed.
- Aluminium-drawer on roller runners $w \times d \times h = 56 \times 60 \times 10 \text{ cm}$ (a maximum of eight drawers is possible).
- Additional length and cross dividers (for drawer).
- Stainless steel drawer (liquid tight) $w \times d \times h = 56 \times 60 \times 10 \text{ cm}$ (a maximum of eight drawers is
- Additional shelf on rails or overlays.
- Basket on rails, w x d x h = 60 x 65 x 10 cm.
- Aluminium tray, w x d x h = $60 \times 65 \times 2.5$ cm, with rails.
- Vario Support, Storage Shelf on roller runners.
- Condensate container can be emptied manually if no condensation is required (e.g., in an operating theatre).
- GSM-Module for alarm text message transmission.
- Temperature documentation.
 - PC-KIT-NET (network capable version)
 - 4-pin PT 100 or PT 1000 temperature sensor, class 1/3 B, including built-in reference body (measurement bottle or refrigerating block).
 - Mechanical pen-recording thermometer with two waxed paper strips.