

LABSOLUTE® PIPETTE TIPS COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| LABSOLUTE® Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|------------------------|-----------------------|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Single-channel pipette | 0.1–2 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Single-channel pipette | 2–20 µl | 1 | | | | x | | x | x | x | | | | | | | | | |
| Single-channel pipette | 10–100 µl | 1 | | | | | x | x | x | x | x | | | | | | | | |
| Single-channel pipette | 20–200 µl | 1 | | | | | | x | x | x | x | | | | | | | | |
| Single-channel pipette | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Single-channel pipette | 500–5000 µl | 1 | | | | | | | | | | | | | x | | | | |
| Single-channel pipette | 1–10 ml | 1 | | | | | | | | | | | | | | | | x | |
| Multi-channel pipette | 5–50 µl | 8/12 | | | | | x | | x | x | | | | | | | | | |
| Multi-channel pipette | 20–200 µl | 8/12 | | | | | | x | x | x | x | | | | | | | | |
| Multi-channel pipette | 50–300 µl | 8/12 | | | | | | | | | | x | x | | | | | | |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Sartorius/BioHit Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|-----------------------------|--|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| mLine | 0.1–3 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| mLine | 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| mLine | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| mLine | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | |
| mLine | 5–100 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| mLine | 30–300 µl | 8/12 | | | | | | | | | | x | x | | | | | | |
| Proline | 0.1–2.5 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Proline | 2–20 µl 5–50 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Proline | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Proline | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | |
| Proline | 5–50 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| Proline | 50–300 µl | 8/12 | | | | | | | | | | x | x | | | | | | |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Brand Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|-----------------|--|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Transferpette | 0.1–1 µl 0.5–10 µl | 1 | | | | | | | | | | | | | | | | | |
| Transferpette | 2–20 µl | 1 | | | x | | | | | | | | | | | | | | |
| Transferpette | 2–20 µl 5–50 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Transferpette | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Transferpette | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | |
| Transferpette | 2–20 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| Transferpette | 2.5–25 µl 5–50 µl 10–100 µl 20–200 µl | 8/12 | | | | x | x | x | | x | x | | | | | | | | |
| Transferpette | 30–300µl | 8/12 | | | | | | | | | | x | x | | | | | | |
| Transferpette | 500–5000 µl | 1 | | | | | | | | | | | | | | | x | | |
| Transferpette S | 0.1–1 µl 0.1–2.5 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Transferpette S | 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Transferpette S | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Transferpette S | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | |
| Transferpette S | 5–50 µl 10–100 µl 20–200 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| Transferpette S | 30–300 µl | 8/12 | x | x | | | | | | | | x | x | | | | | | |
| Transferpette S | 500–5000 µl | 1 | | | | | | | | | | | | | | | x | | |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Eppendorf Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|-------------------|-----------------------------------|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Research Plus | 0.1–2.5 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Research Plus | 2–20 µl | 1 | | | x | | | | | | | | | | | | | | |
| Research Plus | 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Research Plus | 30–300 µl | 1 | | | | | | | | | x | x | | | | | | | |
| Research Plus | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Research Plus | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | |
| Research Plus | 10–100 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| Research Plus | 30–300 µl | 8/12 | | | | | | | | | | x | x | | | | | | |
| Research Plus | 500–5000 | 1 | | | | | | | | | | | | | | x | | | |
| Research Plus | 1000–10000 µl | 1 | | | | | | | | | | | | | | | | x | x |
| Reference | 0.1–2.5 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Reference | 2–20µl | 1 | | | x | | | | | | | | | | | | | | |
| Reference | 2–20 µl 10–100 µl 50–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Reference | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Reference | 500–5000 µl | 1 | | | | | | | | | | | | | x | | | | |
| Reference | 1000–10000 µl | 1 | | | | | | | | | | | | | | | | x | x |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Thermo/Finn Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|------------------------|--|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Finnpette 4500/ F1/F2 | 0.2 – 2 µl 0.5 – 5 µl | 1 | x | | x | | | x | x | x | | x | | | | | | | |
| Finnpette 4500/ F1/F2 | 0.5 – 5 µl 2 – 20 µl 5 – 50 µl 10 – 100 µl 20 – 200 µl | 1 | | | | x | x | x | x | x | x | x | | | | | | | |
| Finnpette 4500/ F1/F2 | 100 – 1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Finnpette 4500/ F1/F2 | 0.5 – 10 µl | 8/12 | x | | x | | | x | x | x | | x | | | | | | | |
| Finnpette 4500/ F1/F2 | 5 – 50 µl | 8/12 | | | | x | x | x | x | x | x | x | | | | | | | |
| Finnpette 4500/ F1/F2 | 50 – 300 µl | 8/12 | | | | x | x | x | x | x | x | x | | | | | | | |
| Finnpette 4500/ F1/F2 | 500 – 5000 µl | 1 | | | | | | | | | | | | | | | x | | |
| Finnpette 4500/ F1/F2 | 1000 – 10000 µl | 1 | | | | | | | | | | | | | | | | | x |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Gilson Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|----------------------------------|-----------------------------------|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Pipetman Classic, G, L: P2 | 0.2–2 µl 1–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Pipetman Classic, G, L: P20 | 2–20 µl 20–100 µl 50–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Pipetman Classic, G, L: P1000 | 200–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Pipetman P10L | 0.5–10 µl 2–20 µl 20–200 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| Pipetman P300L | 20–300 µl | 8/12 | | | | | | | | | | x | x | | | | | | |
| Pipetman Classic, L, G | 500–5000 µl | 1 | | | | | | | | | | | | | | x | | | |
| Pipetman Classic, L, G | 1000–10000 µl | 1 | | | | | | | | | | | | | | | | | x |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Nichiryo Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|----------------------------------|-----------------------------------|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Nichipet Premium, EX, EX Plus | 0.1–2 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Nichipet Premium, EX, EX Plus | 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Nichipet Premium, EX, EX Plus | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Nichipet 7000 | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | |
| Nichipet 7000 | 5–50 µl 40–200 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | |
| Nichipet 7000 | 50–300 µl | 8/12 | | | | | | | | | | x | x | | | | | | |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Rainin Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | |
|---|-----------------------------------|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 |
| Pipet-Lite SL-2XLS | 0.1–2 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Pipet-Lite SL-20XLS | 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Pipet-Lite SL-300XLS | 20–300 µl | 1 | | | | | | | | | | x | x | | | | | | |
| Pipet-Lite SL-1000XLS Classic PR-2 Classic PR-20 Classic PR-1000 | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Classic PR-2 | 0.1–2 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | |
| Classic PR-20 | 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | |
| Classic PR-1000 | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | |
| Classic PR-5000 | 500–5000 µl | 1 | | | | | | | | | | | | | | x | | | |
| Classic PR-10000 | 5000–10000 µl | 1 | | | | | | | | | | | | | | | | | x |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021

LABSOLUTE® PIPETTE TIPS
COMPATIBILITY CHART



We confirm that LABSOLUTE® pipette tips fit on following pipette models*:

| Socorex Pipette | Volume range | Compatibility: Single-channel (1), multi-channel (8/12) pipette | 10 µl | | 10 µl XL/ 20 µl | 20 µl | 100 µl | 200 µl | | | 300 µl | | 1,000 µl | 1,250 µl | 5 ml | | 10 ml | | | |
|-----------------|---|--|---------------------------------------|---------|---|---------|---------|------------------|------------------|---|---------|---------------------------------------|----------|---|--|---------|---------|------------------|---------|---------|
| | | | 7695840, 7695860, 7695870, 7695890 | 7695880 | 7695841, 7695861, 7695871, 7695891, 7695881 | 7695882 | 7695883 | 7695842, 7695843 | 7695844, 7695845 | 7695862, 7695863, 7695872, 7695873, 7695892 | 7695884 | 7695846, 7695864, 7695874, 7695893 | 7695885 | 7695847, 7695848, 7695865, 7695866, 7695875, 7695876, 7695894, 7695886 | 7695849, 7695867, 7695877, 7695878, 7695887, 7695854, 7695868 | 7695850 | 7695933 | 7695855, 7695856 | 7695851 | 7695934 |
| Calibra® 822 | 0.2–2 µl 1–10 µl | 1 | x | x | x | | | | | | | | | | | | | | | |
| Calibra® 822 | 1–10 µl 2–20 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | | |
| Calibra® 822 | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | | |
| Calibra® 852 | 1–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | | |
| Calibra® 852 | 10–100 µl 20–200 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | | |
| Acura 825 | 0.1–2 µl 0.5–10 µl | 1 | x | x | x | | | | | | | | | | | | | | | |
| Acura 825 | 1–10 µl 2–20 µl 5–50 µl 10–100 µl 20–200 µl | 1 | | | | x | x | x | x | x | x | | | | | | | | | |
| Acura 825 | 100–1000 µl | 1 | | | | | | | | | | | x | x | | | | | | |
| Acura 855 | 0.5–10 µl | 8/12 | x | x | x | | | | | | | | | | | | | | | |
| Acura 855 | 5–50 µl 20–200 µl | 8/12 | | | | x | x | x | x | x | x | | | | | | | | | |
| Acura 855 | 40–350 µl | 8/12 | | | | | | | | | | | x | | | | | | | |
| Acura 835 | 500–5000 µl | 1 | | | | | | | | | | | | | | | | | x | |
| Acura 835 | 1000-10000 µl | 1 | | | | | | | | | | | | | | | | | | x |

* There is no claim on completeness, actuality, quality and accuracy. Status: August 2021