


InLab[®] 741-ISM

Cond °C ISM



For low conductivity range samples like pure water and ultra pure water.



Sensor type conductivity cell

Order number 30014094

Specifications

measuring range 0.001...500 µS/cm

zero point (pH) -

temperature range 0...100 °C

temperature probe NTC 30 kΩ

shaft material stainless steel

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.2 m

connector Mini-DIN

cell constant (cm⁻¹) 0.105 (certified)

cell type 2 steel poles

indicating metal -

Storage

dry