


# InLab<sup>®</sup> Redox Ag850

ORP half-cell with silver ring

Scale 1:1



**Sensor type** ORP half-cell

**Order number** 59904408

**Specifications**

pH range -  
zero point (pH) -

temperature range -30...80 °C  
temperature probe -

shaft material polypropylene  
membrane glass -  
membrane resistance -

reference system -  
reference diaphragm -  
reference electrolyte -  
bridge electrolyte -

electrode head S7  
cable length -  
connector -

cell constant (cm<sup>-1</sup>) -  
cell type -

indicating metal silver tip

**Storage**

dry