


InLab[®] Mono

pH half-cell

pH Ref °C ISM



Scale 1:1



Sensor type pH half-cell

Order number 51343195

Specifications

pH range 0...14

zero point (pH) 7.00 ± 0.25

temperature range 0...100 °C

temperature probe -

shaft material glass

membrane glass HA

membrane resistance < 600 MΩ

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)