


InLab® Solids Pro-ISM

pH Ref °C ISM



For solids and semisolids samples like meat, cheese, fruits, vegetables

Scale 1:1



Sensor type combined pH electrode

Order number 51344155

Specifications

pH range 1...11

zero point (pH) 7.00 ± 0.25

temperature range 0...80 °C

temperature probe NTC 30 k Ω

shaft material glass

membrane glass Lot

membrane resistance < 250 M Ω

reference system ARGENTHAL™

reference diaphragm open junction

reference electrolyte XEROLYT® EXTRA polymer

bridge electrolyte -

electrode head MultiPin™

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)

or FRISCOLYT-B® (51340053)