

# Accessories

| Accessories for Universal Platform, 46 × 46 cm, S/S (Item No. 30400054) | Description                           | Max per Platform |
|---|---------------------------------------|------------------|
| 30400084  | Flask Clamp, 10 mL                    | 113              |
| 30400085  | Flask Clamp, 25 mL                    | 64               |
| 30400086  | Flask Clamp, 50 mL                    | 32               |
| 30400087  | Flask Clamp, 125 mL                   | 20               |
| 30400088  | Flask Clamp, 250 mL                   | 20               |
| 30400089  | Flask Clamp, 500 mL                   | 13               |
| 30400090  | 1 Liter Erlenmeyer Flask Clamp        | 8                |
| 30400091  | 2 Liter Erlenmeyer Flask Clamp        | 5                |
| 30400092  | 2.8 Liter Erlenmeyer Flask Clamp      | 2                |
| 30400093  | 3 Liter Erlenmeyer Flask Clamp        | 2                |
| 30400094  | 4 Liter Erlenmeyer Flask Clamp        | 4                |
| 30400095  | 5 Liter Erlenmeyer Flask Clamp        | 4                |
| 30400096  | 6 Liter Erlenmeyer Flask Clamp        | 2                |
| 30400097  | 500 mL Media Bottle Clamp             | 16               |
| 30400098  | 1000 mL Media Bottle Clamp            | 10               |
| 30400099  | Flask Clamp PVC 125 mL                | 20               |
| 30400100  | Flask Clamp PVC 250 mL                | 18               |
| 30400101  | Flask Clamp PVC 500 mL                | 12               |
| 30400102  | Flask Clamp PVC 1000 mL               | 8                |
| 30400103  | Flask Clamp PVC 2000 mL               | 4                |
| 30400105  | Test Tube Rack 13 mm Pivoting         | 2                |
| 30400106  | Test Tube Rack 16 mm Pivoting         | 2                |
| 30400107  | Test Tube Rack 20 mm Pivoting         | 2                |
| 30400108  | Test Tube Rack 25 mm Pivoting         | 2                |
| 30400109  | Test Tube Rack 30 mm Pivoting         | 2                |
| 30400110  | Test Tube Rack 10-14 mm Stationary    | 3                |
| 30400111  | Test Tube Rack 16-20 mm Stationary    | 3                |
| 30400112  | Test Tube Rack 21-25 mm Stationary    | 3                |
| 30400113  | Test Tube Rack 50 mL Tubes Stationary | 3                |
| 30400114  | Micro Tube Rack                       | 4                |
| 30400115  | Test Tube Rack 10-13 mm Diameter      | 4                |
| 30400116  | Test Tube Rack 14-16 mm Diameter      | 4                |
| 30400117  | Test Tube Rack 18-20 mm Diameter      | 4                |
| 30400118  | Test Tube Rack 22-25 mm Diameter      | 4                |
| 30400119  | Test Tube Rack For 15 mL Tubes        | 4                |
| 30400120  | Test Tube Rack For 50 mL Tubes        | 4                |
| 30400104  | Clamp Microplate Stainless Steel      | 12               |
| 30400062  | Rubber Mat, 46 × 46 cm                | 1                |



30400054  
Universal Platform, 46 × 46 cm



30400088  
Flask Clamp, 250 mL



30400096  
6 Liter Erlenmeyer Flask Clamp