

Accessories



3011




3012



3013




3014



Order No. 3965

Non-slip mat
for the shaking platform,
420 x 420 mm size, for slow
moving, e.g. nutrient solutions
in Petri dishes



Order No. 3967

Universal mount
for secure attachment of
different shaking objects
between the six rubber-
coated bars



Order No. 3966

Shaking tray
made of stainless steel,
450 x 450 mm, with holes to
accommodate clamps for
Erlenmeyer flasks and other
accessory equipment



Order No. 3968

Platform frame
with four levels, with top rims,
three levels removable. Frame
and trays made of stainless steel.
Distance between trays - 60 mm.
Grooved rubber mats on trays to
ensure secure grip of objects.



Clamps
made of stainless steel for
Erlenmeyer flasks, for screwing
on to shaking tray 3966,
complete with attachment
materials

- Order No. 3983 for 25 ml flasks (79*)
- Order No. 3984 for 50 ml flasks (49*)
- Order No. 3985 for 100 ml flasks (36*)
- Order No. 3986 for 200 ml flasks (22*)
- Order No. 3987 for 250-300 ml flasks (16*)
- Order No. 3988 for 500 ml flasks (12*)
- Order No. 3989 for 1000 ml flasks (9*)
- Order No. 3990 for 2000 ml flasks (4*)

Other clamps on request

* Max. number of clamps per shaking tray