

Accessories



3015



3017



3016



3018



Shaking tray
made of stainless steel,
450 x 450 mm, with holes
to accommodate clamps for
Erlenmeyer flasks and other
accessory equipment

Order No. 3966




Clamps
made of stainless steel for
Erlenmeyer flasks, for screwing
on to shaking tray 3966,
complete with attachment
materials

- Order No. 3983 for 25 ml flasks (79*)
- Order No. 3984 for 50 ml flasks (49*)
- Order No. 3985 for 100 ml flasks (36*)
- Order No. 3986 for 200 ml flasks (22*)
- Order No. 3987 for 250-300 ml flasks (16*)
- Order No. 3988 for 500 ml flasks (12*)
- Order No. 3989 for 1000 ml flasks (9*)
- Order No. 3990 for 2000 ml flasks (4*)

Other clamps on request

* Max. number of clamps per shaking tray



Test tube rack
made of stainless steel, for test
tubes; perforated section can be
tilted by 90°, with clamping
springs for secure grip and silent
shaking of glasses; can be
screwed on to shaking tray 3966

Order No. 3953
for a max. of 24 tubes: 12 - 17 mm diameter
capacity of tray: six test tube racks




Clamps
made of stainless steel, for
separating funnels, for screwing
on to shaking tray 3966

- Order No. 3957 for 50 ml separating funnels (11*)
(ISO/Squibb, neck 19/26)
- Order No. 3958 for 100 ml separating funnels (11*)
(ISO/Squibb, neck 19/26)
- Order No. 3959 for 250 ml separating funnels (8*)
(ISO/Squibb, neck 19/26)
- Order No. 3955 for 250 ml separating funnels (8*)
(conical form, neck 29/32)
- Order No. 3956 for 500 ml separating funnels (6*)
(ISO/Squibb, neck 29/32)

Other clamps on request
* Max. number of clamps per shaking tray

Holding Frame
made of stainless steel, for attachment of a test plate, for screwing
on to shaking tray 3966, complete with attachment materials
pict. see Shakers 3019 - 3020

Order No. 3910
Capacity of tray: six test plates



Non-slip mat
for the shaking platform,
420 x 420 mm size, for slow
moving, e.g. nutrient solutions
in Petri dishes

Order No. 3965



Universal mount
for secure attachment of
different shaking objects
between the six rubber-
coated bars

Order No. 3967