



# Accessories



3019



3020

## Shaking tray

made of 3-mm-thick anodised aluminium, 670 x 537 mm, with holes to accommodate clamps for Erlenmeyer flasks and other accessory equipment



Order No. 3980

## Clamps

made of stainless steel for Erlenmeyer flasks, for screwing on to shaking tray 3980, complete with attachment materials




- Order No. 3983 for 25 ml flasks (99\*)
- Order No. 3984 for 50 ml flasks (99\*)
- Order No. 3985 for 100 ml flasks (50\*)
- Order No. 3986 for 200 ml flasks (26\*)
- Order No. 3987 for 250-300 ml flasks (26\*)
- Order No. 3988 for 500 ml flasks (26\*)
- Order No. 3989 for 1000 ml flasks (12\*)
- Order No. 3990 for 2000 ml flasks ( 9\*)

Other clamps on request

\* Max. number of clamps per shaking tray

## Test tube rack

made of stainless steel for test tubes; perforated section can be tilted by 90°, with clamping springs for secure grip and silent shaking of glasses; can be screwed on to shaking tray 3980



Order No. 3953

for a max. of 24 tubes: 12 - 17 mm diameter  
Capacity of tray: nine test tube racks



Order No. 3910

Capacity of tray: fifteen test plates

## Holding frame

made of stainless steel, for attachment of a test plate, for screwing on to shaking tray 3980, complete with attachment materials



Order No. 3982

## Universal mount

for secure attachment of shaking objects. The mount consists of two longitudinal bars and six transverse bars. Screw-type clamps are used for adjusting the longitudinal bars vertically and transverse bars horizontally.



Order No. 3981

## Rack frame

for doubling the capacity by accommodating two shaking trays 3980; distance between trays - 205 mm. The bottom tray can be loaded with flasks up to a max. of 500 ml.