




Digital conductivity sensor ideal for low conductivity range samples like pure water and ultra pure water.



Sensor type conductivity cell

Order number 30014097

Specifications

measuring range 0.0001...1000 µS/cm

zero point (pH) -

temperature range 0...100 °C

temperature probe Pt1000

shaft material titanium

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.8 m

connector Mini-LTW

cell constant (cm⁻¹) 0.01 (certified)

cell type 2 titanium poles

indicating metal -

Storage

dry