

# Accessories for Immersion circulators IKA ICC



CS.ICC

CL.ICC



BS.ICC

PCS.ICC

Product	Description	Ident. No.
<b>ACCESSORIES FOR ICC</b>		
CS.ICC	Cover, small for ICC, 8 / 9 l bath vessels	0004471500
CM.ICC	Cover, medium for ICC, 12 l bath vessels	0025000290
CL.ICC	Cover, big for ICC, 18 / 20 l bath vessels	0004471600
CLG.ICC	Gable lid for ICC, bath size L	0020023680
CLR.ICC	Cover with 6 ring openings for ICC, 18/20 l bath vessels	0025006276
BS.ICC	Bridge, small for ICC, 8 / 9 l bath vessels	0020003077
BL.ICC	Bridge, big for ICC, 18 / 20 l bath vessels	0020003078
BS.ICC lite	Bridge, small for ICC lite, 8/9 l bath vessels	0020119539
BL.ICC lite	Bridge, big for ICC lite, 8/9 l bath vessels	0020119540
CC2	Cooling coil for ICC	0025001061
PCS.ICC	Pump connection set for ICC	0004471900

## BATH VESSELS

### Bath size S


IB 8 eco	Bath vessel, plastic, 8 l	0004248100
IB R 9 eco	Bath vessel, plastic, 9 l, rectangular	0020004382
IB 9 pro	Bath vessel, stainless steel, 9 l	0004248500

### Bath size M


IB 12 pro	Bath vessel, stainless steel, 12 l	0004577500
-----------	------------------------------------	------------

### Bath size L

IB 18 eco	Bath vessel, plastic, 18 l	0004248200
IB HF 24 eco	Bath vessel, plastic, 24 l	0020006884
IB R 20 eco	Bath vessel, 20 l, plastic, rectangular	0020004383
IB R 52 eco	Bath vessel, 52 l, plastic, rectangular	0020007044
IB R RO 15 eco	Bath vessel for magn. stirrer RO 15	0020006341
IB R RO 15 pro	Bath vessel for magn. stirrer RO 15	0025007761
IB 20 pro	Bath vessel, stainless steel, 20 l	0004248600



IB R 9 eco



IB 12 pro



IB HF 24 eco



IB 20 pro



HF.Si.20.200.50



LT 5.32



H.SI.8



H.FKM.12.ISO



ISO.8



ISO.12

Product	Description	Ident. No.
<b>THERMOFLUIDS</b>		
HF.Si.20.250.50A	Heating fluid Si250A; 10 kg	0020003521
HF.Si.20.200.50	Heating fluid Si200; 10 kg	0020003520
UF.Si.N30.150.10LV	Universal fluid Si150 LV; 9 kg	0020003518

## HOSES

LT 5.32	Temperature control hose, 1.6 m	0020113271
H.PVC.8	Hose, PVC, Ø 8 mm, 2 x 1.5 m	0004568800
H.PVC.12	Hose, PVC, Ø 12 mm, 2 x 1.5 m	0004568900
H.SI.6	Hose, silicone, Ø 6 mm, 3 m	0020102324
H.SI.8	Hose, silicone, Ø 8 mm, 2 x 1.5 m	0004569000
H.SI.12	Hose, silicone, Ø 12 mm, 2 x 1.5 m	0004569100
H.PUR.8.R	Tube PUR 90 °C, Ø 8, 2 x 1.5 m, reinforced	0020004612
H.PUR.12	Tube PUR 90 °C, Ø 12, 2 x 1.5 m	0020004613
H.FKM.8	Tube FKM, Ø 8 mm, -30 – 180 °C, 2 x 1.5 m	0004569200
H.FKM.12	Tube FKM, Ø 12 mm, -30 – 180 °C, 2 x 1.5 m, M 16 x 1	0004569300
H.FKM.8.ISO	Hose set FKM incl. insulation, Ø 8 mm	0010011653
H.FKM.12.ISO	Hose set incl. insulation, FKM, Ø 12 mm	0010011652
H.PUR.12.AD	Hose set incl. M 16x1 connection, PUR 90 °C, Ø 12 mm	0010011654

## INSULATION FOR HOSES

ISO.8	Insulation, Ø 8 mm, 2 x 1.5 m	0004569400
ISO.12	Insulation, Ø 12 mm, 2 x 1.5 m	0004569500


# Accessories for Immersion circulators IKA ICC



Tube rack, 13 mm, S, stainless



Variable rack, ICC, ML, stainless



Variable bottom, ICC, L, stainless



Bottom plate L



Floating tube rack 1



Floating globes, PP



WH10

Product	Description	Ident. No.
<b>TUBE RACKS</b>		
Tube rack, 13 mm, S, stainless	Stainless steel rack for bath size S, Ø 13 mm	0020004026
Tube rack, 17 mm, S, stainless	Stainless steel rack for bath size S, Ø 17 mm	0020004027
Tube rack, 22 mm, S, stainless	Stainless steel rack for bath size S, Ø 22 mm	0020004028
Tube rack, 13 mm, ML, stainless	Stainless steel rack for bath size M+L, Ø 13 mm	0020004029
Tube rack, 17 mm, ML, stainless	Stainless steel rack for bath size M+L, Ø 17 mm	0020004030
Tube rack, 22 mm, ML, stainless	Stainless steel rack for bath size M+L, Ø 22 mm	0020004031
Variable rack, ICC, ML, stainless	Variable rack for bath size M+L.ICC	0020004032
Inlay 1, variable rack, ICC, ML	Inlay for variable rack, Ø 13 mm, 2 pcs.	0020004033
Inlay 2, variable rack, ICC, ML	Inlay for variable rack, Ø 17 mm, 2 pcs.	0020004034
Inlay 3, variable rack, ICC, ML	Inlay for variable rack, Ø 22 mm, 2 pcs.	0020004035
Variable bottom, ICC, L, stainless	Variable bottom for bath size L.ICC	0020004614
Bottom plate L	Bottom plate for bath size L	0020006212
Bottom plate L, ICC	Bottom plate for bath size L.ICC	0020007353
<b>FLOATING RACKS</b>		
Floating tube rack 1	Floating tube rack for 24 x 1.5 / 2.0 ml tubes, 5 pcs.	0020003667
Floating tube rack 2	Floating tube rack for 8 x 15 ml tubes, 5 pcs.	0020003668
Floating tube rack 3	Floating tube rack for 4 x 50 ml tubes, 5 pcs.	0020003669
<b>ADDITIONAL ACCESSORIES</b>		
WH10	WiCo Holder	0020000984
PT 100.30	Temp. sensor, stn. steel	0004284700
PT 100 extension (Lemo)	Extension cable PT 100, lemo connector, 3 m	0020004629
Floating globes, PP	Floating globes, PP, 500 pcs.	0020003666
Distribution Pipe	4-sided distributor tube for improved flow in large tanks, stainless steel	0020027440
PDS.IB.R.RO	Pump connection set ICC for IB R RO baths, incl. distribution pipe	0020026449



CO V 1



Ball valve M16x1

Mechanic fluid level controller



Elbow tube 90°



Olive for DN 6 hoses



Adapter NTP 1/4



Adapter NTP 3/4

Product	Description	Ident. No.
<b>VALVES AND LEVEL CONTROLLERS</b>		
MV 1	Magnetic valve for cooling water control, max. 100 °C	0020003763
CO V 1	Closed pressure valve for external pump connection set, max. 180 °C	0020000249
Ball valve M16 x 1	Ball valve, completely opening (M 16 x 1)	0020004620
Mechanic fluid level controller	Mechanic fluid level controller	0020004618
<b>ADAPTERS</b>		
Elbow tube 90°	Elbow tube 90°	0025001212
Elbow Olive 90° for pump	90° tube adapter for pump mounting	0025007462
Y-hose connector 8 mm	Splitter, PP, max. 135 °C	0020010963
Olive for DN 6 hoses	Olive for DN 6 hoses, 2 pcs.	0020004667
Olive for DN 8 hoses	Olive for DN 8 hoses, 2 pcs.	0020004566
Olive for DN 10 hoses	Olive for DN 10 hoses, 2 pcs.	0020004568
Olive for DN 12 hoses	Olive for DN 12 hoses, 2 pcs.	0020004889
Adapter NTP 1/4	Adapter M 16 to NPT 1/4, 2 pcs.	0020004569
Adapter NTP 1/2	Adapter M 16 to NPT 1/2, 2 pcs.	0020004570
Adapter NTP 3/4	Adapter M 16 to NPT 3/4, 2 pcs.	0020004571
Adapter M16 x 1 f - G1/4 m	Adapter M 16 x 1 female to G 1/4 male	0020008807
Adapter M16 x 1 m - G1/4 m	Adapter M 16 x 1 male to G 1/4 male	0020008734
Adapter M16 x 1 m - G1/2 m	Adapter M 16 x 1 male to G 1/2 male	0020008796
Adapter M16 x 1 f - G1/2 m	Adapter M 16 x 1 female to G 1/2 male	0020008800
Adapter M16 x 1 f - G3/8 m	Adapter M 16 x 1 female to G 3/8 male	0020008730
Adapter M16 x 1 m - G3/8 m	Adapter M 16 x 1 male to G 3/8 male	0020008754
Adapter M16 x 1 m - M 16 x 1 m	Adapter M 16 x 1 male to M 16 x 1 male, for hose extension	0025007118
Lock nut	Lock nut, 2 pcs.	0020004583
Stopper	Stopper, 2 pcs.	0020004584