

# ROTARY EVAPORATORS IKA RV | Accessories



RV 10.1

RV 10.11 enlarged surface



RV 10.3

RV 10.4



RV 10.400

RV 10.500



RV 10.600

RV 10.606



RV 10.84

RV 10.301


Product	Description	Ident. No.
<b>ACCESSORIES RV 10, RV 8, RV 3</b>		
RV 10.1	Set of glassware, vertical	0003655300
RV 10.10	Set of glassware, vertical, coated	0003755400
RV 10.11 enlarged surface	Vertical-condenser, enlarged surface	0020019652
RV 10.110 enlarged surface	Vertical-condenser, coated, enlarged surface	0020019654
RV 10.2	Set of glassware, diagonal	0003755300
RV 10.20	Set of glassware, diagonal, coated	0003755500
RV 10.3	Vertical-intensive condenser, manifold	0003741000
RV 10.30	Vertical-intensive condenser, manifold, coated	0003741100
RV 10.4	Dry ice condensor	0003742000
RV 10.40	Dry ice condensor, coated	0003742100
RV 10.5	Vertical-condenser, reflux	0003743000
RV 10.50	Vertical-condenser, reflux, coated	0003743100
RV 10.6	Vertical-intensive condenser, reflux	0003744000
RV 10.60	Vertical-intensive condenser, reflux, coated	0003744100
RV 05.10	Balljoint clamp	0000736600
RV 10.88	Clamp, NS 29	0003683300
RV 10.400	Evaporation cylinder, NS 29/32, 500 ml	0003738800
RV 10.401	Evaporation cylinder, NS 29/32, 1500 ml	0003738900
RV 10.500	Foam brake, NS 29/32	0003739200
RV 10.600	Distilling spider, NS 29/32, 6 sleeves	0003739400
RV 10.601	Distilling spider, NS 29/32, 12 sleeves	0003739500
RV 10.602	Distilling spider, NS 29/32, 20 sleeves	0003739600
RV 10.610	Distilling sleeve, 20 ml	0003919400
RV 10.606	Distilling spider with 5 flasks, 50 ml	0003740800
RV 10.607	Distilling spider with 5 flasks, 100 ml	0003740900
RV 10.80	Evaporation flask, NS 29/32, 50 ml	0003740100
RV 10.800	Evaporation flask, coated, NS 29/32, 50 ml	0020005330
RV 10.81	Evaporation flask, NS 29/32, 100 ml	0003740200
RV 10.810	Evaporation flask, coated, NS 29/32, 100 ml	0020005331
RV 10.82	Evaporation flask, NS 29/32, 250 ml	0003740300
RV 10.820	Evaporation flask, coated, NS 29/32, 250 ml	0020005332
RV 10.83	Evaporation flask, NS 29/32, 500 ml	0003740400
RV 10.830	Evaporation flask, coated, NS 29/32, 500 ml	0020005333
RV 10.84	Evaporation flask, NS 29/32, 1000 ml	0003740500
RV 10.840	Evaporation flask, coated, NS 29/32, 1000 ml	0004221400
RV 10.85	Evaporation flask, NS 29/32, 2000 ml	0003740600
RV 10.850	Evaporation flask, coated, NS 29/32, 2000 ml	0020005334
RV 10.86	Evaporation flask, NS 29/32, 3000 ml	0003740700
RV 10.860	Evaporation flask, coated, NS 29/32, 3000 ml	0020005335
RV 10.90	Evaporation flask, NS 24/29, 50 ml	0003741200
RV 10.91	Evaporation flask, NS 24/29, 100 ml	0003741300



RV 10.103



RV 10.74



RV 10.8001

RV 10.3000

Product	Description	Ident. No.
<b>ACCESSORIES RV 3, RV 8, RV 10</b>		
RV 10.100	Receiving flask, KS 35/20, 100 ml	0003742200
RV 10.101	Receiving flask, KS 35/20, 250 ml	0003742300
RV 10.102	Receiving flask, KS 35/20, 500 ml	0003742400
RV 10.103	Receiving flask, KS 35/20, 1000 ml	0003742500
RV 10.104	Receiving flask, KS 35/20, 2000 ml	0003742600
RV 10.105	Receiving flask, KS 35/20, 3000 ml	0003742700
RV 10.200	Receiving flask, coated, KS 35/20, 100 ml	0003743200
RV 10.201	Receiving flask, coated, KS 35/20, 250 ml	0003743300
RV 10.202	Receiving flask, coated, KS 35/20, 500 ml	0003743400
RV 10.203	Receiving flask, coated, KS 35/20, 1000 ml	0003743500
RV 10.204	Receiving flask, coated, KS 35/20, 2000 ml	0003743600
RV 10.205	Receiving flask, coated, KS 35/20, 3000 ml	0003743700
RV 10.300	Powder flask, NS 29/32, 500 ml	0003738200
RV 10.301	Powder flask, NS 29/32, 1000 ml	0003738300
RV 10.302	Powder flask, NS 29/32, 2000 ml	0003738400
RV 10.70	Vapor tube, NS 29/32	0003812200
RV 10.74	Vapor tube, short, NS 29/32	0003971700
RV 10.8001	Seal, PTFE, stn. steel	0003907000
RV 10.3000	Extension plate	0003859000

## ACCESSORIES RV 3

RV 3.1	Support pole	0025004587
RV 3.4	Woulff bottle	0025004588
RV 10.5005	Hose kit	0020000552

## ACCESSORIES RV 8

RV 8.1	Support pole	0004438500
RV 8.2	Mechanical end stop for RV 8	0025001342
RV 8.3	Stand pole	0025004707
RV 10.5005	Hose kit	0020000552

## ACCESSORIES VALVE RV 10 AUTO

RV 10.5002	Filter	0003903800
RV 10.5003	Pressure regulating valve	0003907100



RV 3.1

RV 3.4



RV 10.5001



RV 10.5003